



## **MPIC 2.3**

# Release Notes

---

January 2022

## Copyright notice

Copyright © 2022 eMDs, a brand of CompuGroup Medical. All Rights Reserved.

Use of this documentation and related software is governed by a license agreement and this document is incorporated into and is made part of the license agreement. This documentation and related software contain confidential, proprietary, and trade secret information of eMDs, a brand of CompuGroup Medical, and is protected under United States and international copyright and other intellectual property laws. Use, disclosure, reproduction, modification, distribution, or storage in a retrieval system in any form or by any means is prohibited without the prior express written permission of eMDs, a brand of CompuGroup Medical. This documentation and related software is subject to change without notice.

## Publication date

January 2022

## Product

MPIC

## Corporate address

eMDs, Inc.  
10901 Stone Lake Blvd.  
Austin, Texas 78759  
512-257-5200

# Table of contents

<b>Chapter 1 - Enhancements</b> .....	<b>1</b>
SQL Server .....	1
Pending Transactions for CGM LYTEC .....	1



# Chapter 1 - Enhancements

This chapter presents a high-level description of the following enhancements to MPIC.

- If you are upgrading from MPIC 2.1 to MPIC 2.3, there are several installation and configuration steps that are NOT necessary. Consult the MPIC 2.3 User's Guide for more information.
- If you are using CGM MEDISOFT and are currently using MPIC 2.2, you do not need to upgrade to 2.3. Functionality changes affect only CGM LYTEC.

MPIC 2.3 can be downloaded from here:

CGM MEDISOFT	<a href="https://supportcenter.emds.com/richmond/articles/Richmond_VAR_KBA/Medisoft-and-Medisoft-Clinical-Download-Links">https://supportcenter.emds.com/richmond/articles/Richmond_VAR_KBA/Medisoft-and-Medisoft-Clinical-Download-Links</a>
CGM LYTEC	<a href="https://supportcenter.emds.com/richmond/articles/Richmond_VAR_KBA/Lytec-and-Lytec-MD-Download-Links">https://supportcenter.emds.com/richmond/articles/Richmond_VAR_KBA/Lytec-and-Lytec-MD-Download-Links</a>

## SQL Server

SQL Server Express has been removed from the installer for MPIC.

- For CGM LYTEC, MPIC will continue to use the version of SQL Server that CGM LYTEC uses.
- For CGM MEDISOFT, you can download version 2017 of SQL Server Express here if necessary:

<https://www.microsoft.com/en-us/sql-server/sql-server-2017>

---

If MPIC is already installed and you are using SQL Server 2012 Express or Standard or 2017, you do not need to upgrade your version of SQL Server. However, if you are using SQL Server 2008 or older, you must upgrade your version.

---

You can find basic steps to installing SQL Server Express 2017 in the MPIC 2.3 User's Guide.

## Pending Transactions for CGM LYTEC

MPIC can now consume up to 12 Diagnosis codes per charge line (Practice Partner currently only sends up to 6 codes per charge line in both CGM MEDISOFT and CGM LYTEC). Previously, you were limited to 4 codes per charge line for CGM LYTEC.

